


# ADVANCED ASSESSMENT

## Principles Of Oximetry



# ADVANCED ASSESSMENT

## Principles Of Oximetry

### **AUTHORS**

**Rob Theriault AEMCA, RCT(Adv.), CCP(F)**  
Peel Region Base Hospital

### **REVIEWERS/CONTRIBUTORS**

**Donna Smith AEMCA, ACP**  
Hamilton Base Hospital

# How it works

- ❖ A pulse oximeter (SpO<sub>2</sub> monitor) is a non-invasive device which measures the percentage of hemoglobin saturated with oxygen.
- ❖ It emits red and infrared light through vascular tissue and measures the amount of absorbed light on the other side.
- ❖ Hemoglobin changes its shape depending on whether or not it's carrying oxygen molecules.
- ❖ Light absorption also changes depending on whether the hemoglobin is carrying oxygen.
- ❖ This is how a percentage of oxygen “saturation” is derived.


# How it works

- ❖ 98% of oxygen is transported bound to hemoglobin – the rest is transported dissolved in blood plasma
- ❖ hemoglobin is part of the red blood cell
- ❖ each hemoglobin can carry 4 molecules of O<sub>2</sub> (said to be fully saturated)
- ❖ normal saturation is 95-98% (or as high as 100%) on room air
- ❖ a saturation of 94% = hypoxemia (strict definition of hypoxemia is: PaO<sub>2</sub> less than 80 mmHg)

# How it works

- ❖ tissue displacement is also incorporated into the SpO<sub>2</sub> calculation for accuracy.
- ❖ Tissue displacement is represented on the SpO<sub>2</sub> monitor in the form of a plethysmograph - sometimes called the “pleth” for short or a “pulse waveform”.
- ❖ Note: SpO<sub>2</sub> monitors are sometimes sensitive to motion artifact. If you look at the plethysmograph and the waveform morphology (shape) is the same with each pulse, then the numeric value (% saturation) will be accurate.


# plethysmograph



Note that the shape of the “Pleth” is consistent and that each wave corresponds to an ECG wave above. This tells you that the numeric value is accurate

# How it works - adult


- ❖ In an adult, you generally apply the SpO<sub>2</sub> probe to the finger (see below).
- ❖ Be sure that the light source is directly over the sensor - i.e. be careful not to push the probe on too far. The light source should be directly over the nailbed




Finger probe

# How it works - pediatric


- ❖ In an infant you might choose a disposable probe to place around the big toe.



Disposable probe



Wrap and secure toe so that the light source is directly over the sensor



# Indications for SpO<sub>2</sub> monitoring

- ❖ any patient who requires supplemental O<sub>2</sub>
- ❖ any patient receiving analgesia and/or sedation
- ❖ as an adjunct assessment tool for anyone requiring positive pressure ventilation (PPV) by mask or by ETT
- ❖ to assist in ruling out the need for O<sub>2</sub> therapy
- ❖ any patient to be transported by air

# Benefits of SpO<sub>2</sub> monitoring (cautious)

- ❖ continuous monitoring of oxygenation status
- ❖ to monitor the efficacy of therapeutic interventions. e.g. O<sub>2</sub> therapy, treatment with Salbutamol, NTG for A.P.E., etc
- ❖ guide for efficacy of ventilation-oxygenation of the patient who is receiving PPV
- ❖ guide for oxygenation needs in the patient who is being intubated



# Clinical points

- ❖ Whenever possible, the Paramedic should attempt to obtain an SpO<sub>2</sub> reading on room air, followed by a reading with supplemental oxygen
- ❖ Important: NEVER withhold oxygen from the patient in severe respiratory distress

*A drop in SpO<sub>2</sub> may occur before the development of cyanosis in a patient whose respiratory status is deteriorating – therefore SpO<sub>2</sub> may provide an early warning sign*

# Limitations of of SpO<sub>2</sub> monitoring

- ❖ Anemia - the anemic patient may be 100% saturated, but his/her oxygen carrying capacity will be low due to decreased circulating red blood cells (hemoglobin)
- ❖ O<sub>2</sub> therapy may be indicated to fully saturate the existing hemoglobin and add to the amount of dissolved O<sub>2</sub> in the blood plasma

e.g. various types of anemia's - pregnancy

e.g. hypovolemia

# Limitations of SpO<sub>2</sub> monitoring

## EXAMPLES

- ❖ carboxyhemoglobin - the patient who has CO poisoning may have a high SpO<sub>2</sub> reading, however this is because the pulse oximeter cannot distinguish between CO bound to hemoglobin from O<sub>2</sub> bound to hemoglobin
- ❖ hypoperfusion - inadequate pulsations
- ❖ arterial compression or disruption- i.e. from injury or from application of a blood pressure cuff on the same arm that is used to get the SpO<sub>2</sub> reading
- ❖ ambient light
- ❖ nail polish

# Limitations of SpO<sub>2</sub> monitoring

## EXAMPLES

- ❖ motion artifact
- ❖ patients with chronically impaired oxygenation (e.g. emphysemics) will normally saturate low. - if an emphysemic is saturating at 91% but is otherwise in no distress whatsoever, then 91% should not be taken as seriously as one would take a saturation of 91% in a patient with normally healthy lungs.

# SpO<sub>2</sub> summary

- ❖ monitors oxygenation
- ❖ guide to the efficacy of interventions
- ❖ early warning device
- ❖ has its limitations



# Well Done!

Ontario Base Hospital Group  
Self-directed Education Program